



Fixation rings

TC550 for tcpO₂/CO₂ electrodes

The TC550 Fixation rings allow for flexible and PVC-free fixation of transcutaneous electrodes.

The Fixation ring kit

Fixation rings developed and manufactured by Perimed are available in a Fixation ring kit. The Fixation ring kit comprises 50 Fixation rings (TC 550) and one bottle of Contact liquid (TC560). The TC550 Fixations rings are made of polyethylene plastic (PVC-free) to comply with better patient and user environment.

The TC550 Fixation rings have a small base diameter. This makes it possible to attach the rings to fingers, toes or other curvature parts without risking leakage problems. The adhesive of the TC550 Fixation rings is biocompatible according to ISO10993 and designed to firmly attach to skin but still be easy to remove without damaging fragile skin.

One TC560 Contact liquid bottle is sufficient for the use of 50 TC550 Fixation rings, approximately 3-5 drops of TC560 Contact liquid per ring.



Technical Specifications

Fixation ring material:	Polyethylene plastic (PVC-free) Biocompatible adhesive according to ISO10993
Contact liquid:	Propylene glycol and water